



# MELISSA 72 - Mattress overlay

#### A firm and therapeutic mattress overlay up to and including pressure ulcer category 2

The Melissa 12 is a stable mattress overlay designed to reduce the risk of pressure ulcers up to and including category 2. A 6 cm overlay mattress with pressure relieving properties, this mattress is a good complement to your existing homecare mattress. In addition, with the ability to flip and rotate the core, you get long-term use while maintaining quality.

- HIGH RESILIENT COLD FOAM. The Melissa 12 mattress overlay is made from high resilient cold foam which provides a durable and stable mattress.
- WATERPROOF COVER WITH ELASTIC. The mattress overlay comes with a noiseless RX cover. The cover is supplied with elastic bands at the corners to keep the mattress overlay in place.
- REVERSIBLE FOAM CORE.

The foam core of this mattress overlay can be both flipped and rotated without altering its properties. By doing so, the lifespan of your Melissa 12 mattress overlay will be extended.



### MELISSA 12 - Mattress overlay

#### A GREAT COMPLEMENT TO YOUR HEALTHCARE MATTRESS.

Melissa 12 mattress overlay is designed to provide comfort and pressure relief up to and including pressure ulcer category 2. The two-way stretch hygienic cover is designed to reduce the risk of shear and work in conjunction with the foam core to provide excellent pressure relief.

#### Waterproof cover with good properties

Melissa 12 mattress overlay comes with a cover from the RX series. A series of covers made of 2-stretch polyurethane fabric that is both breathable and noiseless.

The cover has a zipper on 3 or 4 sides (RXI-5) to ease the daily care of the mattress.

All covers have a drip shield covering the zipper to facilitate hygiene.

6 cm cold foam50 kg/m<sup>3</sup>

#### STANDARD FOR THE OBJECTIVE EVALUATION OF THE PRESSURE-RELIEVING PROPERTIES OF MATTRESSES.

In order to obtain an objective assessment of the pressure relieving properties of mattresses, RISE has developed a standard where the actual counterpressure of a mattress is measured by means of an indentational body. The measurement method is described in the Swedish standard SS 876 00 13.

Load capacity Melissa 12 mattress overlav has been tested by RISE according to method SS 876 00 13 and is Large indenter 300 N pressure relieving up to wound category 2. Large indenter 600 N Small indenter 40 N

#### The test results are presented below.

Mattress	Large indenter 300 N	Large indenter 600 N	Small indenter 40 N
Melissa 12, mattress overlay	27	41	44

\*RISE (Research Institutes of Sweden), a Swedish independent research institute - www.ri.se

30,0 - 34,9

45,0 - 52,9

45,0 - 52,9

Extra low

< 30,0

< 45,0

< 45,0

## **PRODUCT SPECIFICATIONS**

Materials	6 cm cold foam 50 kg/m <sup>3</sup>	Us
Sizes	200/205/210 x 80/85/90/105/120 x 6 cm	
	(For other sizes contact Järven Health Care.)	Re
Wound category	Up to and including category 2.	
Mattress weight	approx. 5 kg depending on length and width.	M

	User weight	Up to and including 180 kg.	
	Fire safety	SS 87600 01, normal risk (test method EN 597-1 & 2)	
	Recycling	Mattress with cover is fully recyclable. Sort as combustible waste. Disassembly is not required.	
h.	Marking	Year of manufacture, month, day, EAN, care instructions, contact information, standards met and bed location.	



Bromsvägen 3 | SE-891 60 Örnsköldsvik | SWEDEN Tel: +46 (0)660 - 22 28 50 | Fax: +46 (0)660 - 824 10 www.jarven.se | healthcare@jarven.se



High

40,0 - 44,9

61,0 - 68,9

61,0 - 68,9

Medium

35,0 - 39,9

53,0 - 60,9

53,0 - 60,9