



MELISSA FLEX - Overlay

An affordable mattress overlay that provides pressure relief and support

With the Melissa Flex mattress overlay, you get an affordable mattress overlay designed to reduce the risk of developing pressure sores, category 2. With a thickness of 6 cm and pressure-relieving properties, it is a great complement to your existing medical mattress.

Additionally, the ability to flip and rotate the core ensures long-lasting use while maintaining quality.

- VISCOELASTIC MEMORY FOAM. The Melissa Flex mattress overlay features viscoelastic memory foam with pressure-relieving properties.
- WATERPROOF COVER WITH ELASTIC. The mattress overlay comes with a noiseless RX cover. The cover is supplied with elastic bands at the corners to keep the mattress overlay in place.
- REVERSIBLE FOAM CORE.

The foam core of this mattress overlay can be both flipped and rotated without altering its properties. By doing so, the lifespan of your Melisssa 12 mattress overlay will be extended.



MELISSA FLEX - Overlay

A DURABLE AND PRESSURE RELIEVING MATTRESS OVERLAY

The Melissa Flex mattress topper is designed to provide pressure relief and comfort, and can also be used as a overlay and spacer. The two-stretch hygiene cover is designed to work with the mattress foam core to provide good pressure relief. The cover is also designed to reduce the risk of shear.

Waterproof cover with good properties

Melissa Basic mattress overlay comes with a cover from the RX series. A series of covers made of 2-stretch polyurethane fabric that is both breathable and noiseless. Depending on the choice of cover, the cover has zippers on three or four sides. All covers have a drip shield covering the zipper to facilitate hygiene

6 cm viscoelastic memory foam 35 kg/m3

STANDARD FOR THE OBJECTIVE EVALUATION OF THE PRESSURE-RELIEVING PROPERTIES OF MATTRESSES.

In order to obtain an objective assessment of the pressure relieving properties of mattresses, RISE has developed a standard where the actual counterpressure of a mattress is measured by means of an indentational body. The measurement method is described in the Swedish standard SS 876 00 13.

Melissa Flex mattress overlay has been tested by RISE* according to method SS 876 00 13 and is pressure relieving up to wound category 2.

Load capacity	Extra low	Low	Medium	High
Large indenter 300 N	< 30,0	30,0 - 34,9	35,0 - 39,9	40,0 - 44,9
Large indenter 600 N	< 45,0	45,0 - 52,9	53,0 - 60,9	61,0 - 68,9
Small indenter 40 N	< 45,0	45,0 - 52,9	53,0 - 60,9	61,0 - 68,9

The test results are presented below.

Matress	Large indenter 300N	Large indenter 600N	Small indenter 40N
Melissa Flex, bäddmadrass	27	47	46

 $\ensuremath{^*\!\text{RISE}}$ (Research Institutes of Sweden), a Swedish independent research institute - www.ri.se

PRODUKTSPECIFIKATIONER

Material	6 cm viscoelastic memory foam 35 kg/m³
Sizes	200/205/210 x 80/85/90/105/120 x 6 cm (For other sizes contact Järven Health Care.)
Wound category	Up to and including category 2
Matress weight	approx. 4 kg depending on length and width

	User weight	Up to and including 180 kg.
	Fire safety	SS 87600 01, normalrisk (testmetod EN 597-1 & 2)
2.)	Recycling	Mattress with cover is fully recyclable. Sort as combustible waste. Disassembly is not required.
th	Marking	Year of manufacture, month, day, EAN, care instructions, contact information, standards met and bed location.



Bromsvägen 3 | SE-891 60 Örnsköldsvik | SWEDEN Telefon: 0660 - 22 28 50 | Telefax: 0660 - 824 10 www.jarven.se | healthcare@jarven.se