



MOLLIS BLUE BLACK

Sturdy trolley mattress for pressure relief right from the start

When it comes to care work, it is important to ensure the comfort and health of the patient from the very first moment. Mollis Blue Black has been developed to be a sturdy and safe trolley mattress with pressure-relieving properties. Keep the patient's comfort and health in focus from the beginning of the care process with Mollis Blue Black.

- **SELECTABLE UNDERSIDE.** The Mollis Blue Black trolley mattress can be supplied with a glossy/lacquered underside side for increased sliding ability when moving the mattress sideways, and a matt/ uncoated side for higher friction.
- FULLY WELDED FOR GOOD HYGIENE. To ensure good hygiene, the Mollis Blue Black has a welded cover with no seams or zips.
- OPTIONAL: CUSHIONS.

The Mollis Blue Black can be supplied with a matching cushion that has an integrated fastening strap. Alternatively, an anatomically shaped cushion with viscoelastic foam and fastening straps is also available.







MOLLIS BLUE BLACK

FOCUS ON THE PATIENT'S COMFORT AND HEALTH FROM THE VERY BEGINNING OF THE CARE PROCESS.

Often, the trolley mattress is the first contact the patient has with healthcare. It is also common for the stay on the trolley mattress to be for a longer period of time. To reduce the risk of pressure ulcers occurring, it is of great importance that the mattress is adapted for this.

Mollis Blue and Black has been developed to offer pressure relief right from the start and is also comfortable to lie on.

Thoughtful construction

The uniformly shaped mattress has six handles per long side as standard or three fabric handles per long side to facilitate lateral movement. The mattress can also be ordered without handles if this feature is not required. Mollis Blue Black has a core of two different types of foam which makes the mattress relieving and comfortable to lie on even for longer periods of time.

The bottom consists of 5 cm highly elastic cold foam with a 5 cm of flexible and soft viscoelastic foam on top.

5 cm viscoelastic foam 55 kg/m3

5 cm cold foam 35 kg/m3



To obtain an objective assessment of the pressure-relieving properties of mattresses, RISE has developed a standard where the actual counterpressure of a mattress is measured using indentation bodies. The measurement method is described in Swedish standard SS 876 00 13.

Mollis Blue Black has been tested by RISE according to the SS 876 00 13 and is pressurerelieving up to and including category 2.

The test results are presented below.

Load	Extra low	low	Medium	High
Large indenter 300 N	< 30,0	30,0 - 34,9	35,0 - 39,9	40,0 - 44,9
Large indenter 600 N	< 45,0	45,0 - 52,9	53,0 - 60,9	61,0 - 68,9
Small indenter 40 N	< 45,0	45,0 - 52,9	53,0 - 60,9	61,0 - 68,9

Mattress	Large indenter 300 N	Large indenter 600 N	Small indenter 40 N
Mollis Blue Black	36	74	68

PRODUCT SPECIFICATIONS

Material 5 cm viscoelastic foam 55kg/m3 **User weight** Up to and including 150kg.

5 cm cold foam 35 kg/m3 SS 87600 01, normal risk (test method EN 597-1 & 2) **Fire**

Sizes 200 x 60 x 10 cm Recycling Mattress with cover is fully recyclable. Sorted as (For other sizes contact Järven Health Care.) combustible waste (disassembly not required).

Labelling Year of manufacture, month, day, EAN,

care instructions, contact information, standards met and position in bed.



