



MELISSAFLEX

Preventive healthcare mattress with excellent pressure relief and high comfort

Melissa Flex is a healthcare mattress in the Melissa series. It is developed to be part of the daily care of patients. The mattress is designed to be an aid in the treatment of pressure ulcers up to and including category 2. The documented low pressure is relieving for protruding body parts such as shoulders and heels.

- CAREFULLY ASSEMBLED FOAM CORE. Melissa Flex is a 2-layer mattress with carefully selected foams for comfort and pressure relief.
- UP TO AND INCLUDING WOUND CATEGORY 2.
 Melissa Flex is designed to assist in the treatment of patients with pressure ulcers up to and including Wound Category 2.
- THREE WATERPROOF COVERS OF YOUR CHOICE.
 Melissa Flex is supplied as standard with a noiseless
 RX cover without handles. It can also be supplied
 with optional covers either the RXH cover which
 includes handles or the RXI cover which includes
 both handles and an integrated evacuation function.



MELISSA FLEX

DEVELOPED FOR HIGH COMFORT AND VERY GOOD PRESSURE RELIEF.

Waterpfroof cover with good properties

Melissa Flex comes with a cover from the RX series. A series of covers made of 2-stretch polyurethane fabric that is both breathable and noiseless. Depending on the choice of cover, the cover has zippers on three or four sides. All covers have a drip shield covering the zipper to facilitate hygiene.

Melissa Flex provides very good pressure relief

The foam core consists of two layers of carefully selected foams. Together with our RX cover this creates the conditions for very good pressure relief. With a documented low pressure (46 mm/Hg) for protruding body parts such as heels and

shoulders, no separate heel zone is required.

Melissa Flex offers a seat function

The mattress has been fatigue tested according to standard SS 876 00 11 with good results. The mattress offers a integrated seat function with a raised back position in bed.

6 cm viscoelastic foam 50 kg/m³

8 cm high resilient cold foam 50 kg/m³



In order to obtain an objective assessment of the pressure relieving properties of mattresses, RISE has developed a standard where the actual counterpressure of a mattress is measured by means of an indentational body. The measurement method is described in the Swedish standard SS 876 00 13.

Melissa Flex has been tested by RISE* according to method SS 876 00 13 and is pressure relieving up to wound category 2.

The test resu	lts are preser	nted below.
---------------	----------------	-------------

Load capacity	Extra low	Low	Medium	High
Large indenter 300 N	< 30,0	30,0 - 34,9	35,0 - 39,9	40,0 - 44,9
Large indenter 600 N	< 45,0	45,0 - 52,9	53,0 - 60,9	61,0 - 68,9
Small indenter 40 N	< 45,0	45,0 - 52,9	53,0 - 60,9	61,0 - 68,9

Mattress	Large indenter 300 N	Large indenter 600 N	Small indenter 40 N
Melissa Flex	30	41	46

^{*}RISE (Research Institutes of Sweden), a Swedish independent research institute - www.ri.se

PRODUCT SPECIFICATIONS

Materiala 6 cm viscoelastic foam 50 kg/m³ 8 cm highl resilient cold foam 50 kg/m³.

Sizes 200/205/210 x 80/85/90/105/120 x 14 cm

(For other sizes contact Järven Health Care.)

Wound category Up to and including category 2.

Mattress weight approx. 11 kg depending on length and

width

 $\label{eq:Userweight} \textbf{Up to and including 180 kg}.$

Fire safety SS 87600 01, normal risk (test method EN 597-1 & 2)

Recycling Mattress with cover is fully recyclable.

Sort as combustible waste. Disassembly is not required.

Year of manufacture, month, day, EAN, care

instructions, contact information, standards met and

bed location.



Marking